## Worksheet of the student

Registration No.: 12011241

Date of Performance: 29-03-24

Aim: To investigate the intensity of light coming through two crossed Polaroids and to verify the Malus' law.

Apparatus Required: Light source, Power meter, Polaroids (Polorizers).

Observations table:

Maximum Intensity (Polarizer and analyzer axis are parallel) I<sub>0</sub> = 50

S. No.	Analyzer Rotation with respect to polarizer	Ammeter Reading (I <sub>p</sub> )	Experimental relative intensity I <sub>p</sub> /I <sub>0</sub>	Theoretical relative intensity I <sub>th</sub> /I <sub>0</sub>
0	O°	50		1
2	15°	46.65	0.93369	0.93
3	30°	37.50	0.75	0.75
9	45	25.00	0.5	0.5
(5)	60°	12.50	0-25	0.25
0	75°	3.35	0.067600	0.067(0.06
0	90°	0	0	0
8	105°	3.35	0.067(0.04)	0.067(0.06)
9	120°	12.50	0.25	0.25
(b)	135°	2.5	0.5	0.5
(11)	150°	37.50	0.75	0.75
(1)	165°	46.50	0.933	0.93
(3)	180°	50		

Graph: Experiment relative intensity and Theoretical relative intensity Vs Analyzer rotation

Calculations: Experimental value of Relative Intensity =  $\frac{10}{50} = \frac{50}{50} = 0.93$ ,  $\frac{10}{50} = \frac{3.35}{50} = 0.067$ 

The overlical value of velative  $\frac{1}{2}$   $\frac{$ 

- I) I learned about the Malu's law  $I = I_0 \cos^2 \theta$
- 2) I learned about the position of the axis
  Of the analyzer (0) with respect to the axis
  Of polarizer and the polarization Intensity.
- 3) And i learned about the use of polarizers, Light Source, power meter.
- 4) By this i learned the concept of Polaxization of light.

Signature of student K. G. St

Name of student: Kamboji Gianga Nakhw

Roll No: 46

Reg No: 12011241

